

	Metric	CQC Key Question	2018 Report ²	2019 Report ¹	National Aggregate (England & Wales)	National Aspirational Standard	Comparison to other units
2418 operations	Case ascertainment	Well Led	75%	77%	N/A	None	
2334 operations	Risk-adjusted posterior capsule rupture rate	Effective	0.6%	Case ascertainment below 85% threshold	1.1% ³	None	0% 2.1% Case ascertainment below 85% threshold
1368 operations	Risk-adjusted Visual Acuity Loss	Effective	0.3%	Case ascertainment below 85% threshold	0.9% ³	None	0% 1.5% Case ascertainment below 85% threshold



www.nodaudit.org.uk



KEY ONLY

Display of performance boundaries does vary depending on volume of activity. See FAQs for further information.

Produced by HQIP in partnership with 